



Queen's Square[®] Bladder Stimulator

9 Modes - Hand Held - Battery Operated

Selectable



Vibrating Bladder Emptying Device

Use 3xAA Alkaline Batteries
(Not included)

Model PC9-S

This product can be used to initiate urination and maintain bladder emptying. Adjust cord length and Secure stimulator by placing your hand through the attached cord.

1. Place the Stimulator on the lower abdomen (above the pubic bone) before commencing passing urine and apply gentle pressure to activate the vibration. (If you are unable to exert sufficient pressure to activate, then use the on/off button on the back.)

2. Keep the stimulator vibrating while urine is being passed for about one to two (1-2) minutes after you feel you have finished. Please select different mode for most effective result. Press & hold mode button for 2 seconds (while stimulator is on) to activate. Button remains active for 60 seconds.

On/Off button: Press to switch on. Press again to switch off.

Mode Selection button: To activate, press and hold for 2 seconds while stimulator is on. button remains active for 60 seconds.

If you wish further clinical information, please contact: Clinical Nurse Specialist in Uro-Neurology, The National Hospital for Neurology and Neurosurgery, Queen Square, London WC1 3BG


Reference: British Journal of Urology, Volume 80, Issue 2, August 1997. P234-237



Please select the most effective mode for your need using the mode button. Press & hold mode button for 2 seconds to activate (while stimulator is on). Button remains active for 60 seconds.



Vibrating Modes



Malem Medical, 10 Willow Holt,
Lowdham, Nottingham NG14 7EJ UK
Tel: ++ 44 (0) 115 966 4440
Fax: ++ 44 (0) 115 966 4672
E-mail: malem@malem.co.uk



www.malem.co.uk