




MALEM[®]

Musical or Recorded Message attachment for

Potty / Toilet Training

- Supplied with 2 sensors
- Battery operated
- Trainer Sound Box magnetically attaches to suitable potty



MODEL M0 1 1

A Choice of rewards 20 musical tunes or recorded message



- Choice of 20 tunes or record your own message as a reward after urination (up to 35 seconds after start of urination). Music or message will stop after 30 seconds. It can be stopped sooner by pressing the external Start/Stop button.
- Trainer sound box (62mm x 60mm x 27mm) (87 g) with magnetic clip for easy attachment on underside of potty. POTTY NOT SUPPLIED. Tools not required.
- The short sensor (250mm) is for potty use. Stick underneath potty bowl.
- The long sensor (1000mm) is for toilet use (place in position according to gender e.g. back of toilet for male and front of toilet for female so that it will be in contact with passed urine).
- Battery operated 3 x AAA supplied.

Tune list: 1. London Bridge 2. Lullaby 3. Mary had a little lamb 4. Oh, Susanna 5. Old Macdonald had a farm 6. Pop! Goes the weasel 7. Sally, go round 8. Simple Simon 9. Parade of the wooden soldiers 10. Twinkle, twinkle little star 11. Home sweet home 12. Cuckoo waltz 13. The farmer in the dell 14. For Elise 15. Swan Lake 16. Happy Birthday 17. Silent Night 18. Jingle bells 19. The first noel 20. We wish you a merry Christmas

Malem Medical, 10 Willow Holt, Lowdham Nottingham NG14 7EJ, England, UK

Tel: + 44 (0)115 966 4440 Fax: + 44 (0)115 966 4672 E-mail: malem@malem.co.uk Web: www.malem.co.uk

